Baxter Confidential - INTERNAL USE ONLY		
Part Number: 0719007243	Date: 16APR24	Proofread No.: P01
Designer: R.B	Page: 1 of 16	
Colour Reference:	Black	

BiCart®

Baxter Confidential - INTERNAL USE ONLY		
Part Number: 0719007243	Date: 16APR24	Proofread No.: P01
Designer: R.B	Page: 2 of 16	
Colour Reference:	Black	

THIS PAGE IS INTENTIONALLY LEFT BLANK

Baxter Confidential - INTERNAL USE ONLY		
Part Number: 0719007243	Date: 16APR24	Proofread No.: P01
Designer: R.B	Page: 3 of 16	
Colour Reference:	Black	

EN	Package Leaflet for BiCart® product.....	5
FR	Notice d'utilisation des cartouches BiCart®	7
ES	Manual del operador para el producto BiCart®	10
IT	Manuale operativo prodotto BiCart®	12
RU	Руководство по эксплуатации продукта BiCart®	14

Baxter Confidential - INTERNAL USE ONLY		
Part Number: 0719007243	Date: 16APR24	Proofread No.: P01
Designer: R.B	Page: 4 of 16	
Colour Reference:	Black	

EN

Symbols

+40°C

+104°F

Store below
40° C (104° F)

Recyclable

Baxter Confidential - INTERNAL USE ONLY		
Part Number: 0719007243	Date: 16APR24	Proofread No.: P01
Designer: R.B	Page: 5 of 16	
Colour Reference:	Black	

Package Leaflet for BiCart® product Sodium Hydrogen Carbonate Ph. Eur.

EN

DESCRIPTION

The BiCart product is a polypropylene (PP) cartridge containing dry sodium hydrogen carbonate which enables on-line preparation of a liquid sodium hydrogen carbonate solution used in a hydrogen carbonate-based dialysis fluid.

When the BiCart cartridge is attached to the special holder, water is drawn by the dialysis machine through the cartridge, producing a saturated solution of sodium hydrogen carbonate. The dialysis machine mixes the sodium hydrogen carbonate solution with water and the acid concentrate to produce a hydrogen carbonate-based dialysis fluid.

In some rare cases the hydrogen carboate powder may contain discoloured particles or powder being slightly discoloured when exposed to water. These discolourations originate from the raw materials used during the manufacturing process of the sodium hydrogen carbonate powder. Both the powder and the liquid solution obtained from these BiCart products meet the demands of the European Pharmacopeia. The safety and efficacy of the BiCart product are not affected by these discolourations.

INDICATIONS

To be used in a hydrogen carbonate dialysis treatment for patients suffering from acute renal failure, chronic renal failure or acute intoxication with dialysable substances.

Intended use

The BiCart product is intended for preparation of hemodialysis solutions.

The BiCart product is also intended for preparation of hemodiafiltration and hemofiltration substitution fluids on compatible Gambro dialysis systems.

The BiCart product must always be used with a suitable acid concentrate.

SPECIFICATION

Contents

The BiCart cartridge contains dry sodium hydrogen carbonate powder in compliance with European Pharmacopeia.

Capacity

At a machine setting of 140 mmol/l sodium and 34 mmol/l hydrogen carbonate and together with a suitable acid concentrate, as specified below, up to 32 litres of dialysis fluid can be produced per 100 grams of hydrogen carbonate powder.

Requirements on the acid concentrate

The acid concentrate must be designed for use with a pure sodium hydrogen carbonate concentrate. It must also be suitable for a dilution ratio as recommended by the manufacturer of the dialysis machine to obtain the final dialysis fluid composition.

As an example, when the Gambro concentrate 204 is mixed with a dilution ratio of 1:35, this yields the following concentration in mmol/l:

Sodium	103.0	Potassium	2.0
Calcium	1.75	Magnesium	0.5
Chloride	109.5	Acetic acid	3.0

Together with the sodium hydrogen carbonate concentrate from the BiCart cartridge, the final dialysis fluid will have a sodium concentration of 140 mmol/l and a hydrogen carbonate concentration of 34 mmol/l, at a normal machine setting.

Requirements for water

Dialysis machine inlet water: Inlet water quality must comply with local regulations and if no such regulations are available follow ISO 13959.

CAUTIONS

The information in this package leaflet must be carefully read and understood before using the product.

Do not refill the BiCart cartridge.

The concentrate in the BiCart cartridge must be used within 24 hours after being primed.

The operation of the specially adapted dialysis machine is restricted to persons who are trained in dialysis and who are completely familiar with the operating instructions stated in the operator's manual of the corresponding dialysis machine and in this package leaflet.

Do not use if cap is missing.

WARNINGS

DO NOT use the BiCart product if any damage is found or suspected.

DO NOT use any other type of acid concentrate than the type specified under Specification.

DO NOT initiate dialysis without having followed the instructions in the operator's manual of the corresponding dialysis machine regarding any recommended tests to verify the proper composition of the dialysis fluid.

DO NOT perform any unauthorized modifications or alterations to the BiCart product or to the dialysis machine as this may result in a malfunction or have other serious consequences on the safe operation of the dialysis treatment.

CONTRAINDICATIONS

There are no absolute contraindications for hydrogen carbonate dialysis. The BiCart product must be used only on the direction of a physician who has evaluated all the pertinent features of the product in relation to the individual patient.

ADVERSE REACTIONS

Some undesirable effects related to dialysis can occur, for instance hypotension, nausea, vomiting and cramps.

HANDLING INSTRUCTIONS

For best performance keep the BiCart cartridge at 20 to 25° C (68 to 77° F) for 3-4 hours before use.

For best performance gently shake the BiCart cartridge to release any caking of the powder so it is free flowing before use. See operator's manual of the corresponding dialysis machine.

STORAGE AND SHELF LIFE

Storage

Store below 40° C (104° F)

Shelf life

See the expiry date stated on the product label

WARRANTY AND LIMITATIONS OF LIABILITY

a) The manufacturer warrants that the BiCart product has been manufactured in accordance with its specifications and in compliance with Good Manufacturing Practices, other applicable industry standards and regulatory requirements.

If provided with the lot number of the defective product, the manufacturer will, by replacement, remedy manufacturing defects in the product becoming apparent within the stated expiration period.

b) The warranty under paragraph a) above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, expressed or implied, statutory or otherwise, and there are no warranties of merchantability or other warranties, which extended beyond those described in paragraph a) above. The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the product and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury, or expense arising directly or indirectly from the use of the product, whether as a result of any defect therein or otherwise.


c) The manufacturer shall not be liable for any misuse, improper handling, noncompliance with warnings and instructions, damage arising from events after the manufacturer's release of the product, failure or omission to inspect the product before use in order to insure that the product is in proper condition, or any warranty given by independent distributors or dealers.

This leaflet was prepared by Vantive ULC.
Last revised: 12 April 2024

Baxter Confidential - INTERNAL USE ONLY		
Part Number: 0719007243	Date: 16APR24	Proofread No.: P01
Designer: R.B	Page: 6 of 16	
Colour Reference:	Black	

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Vantive, Bicap and Gampro are trademarks of Vantive Health Inc. or its subsidiaries.
Vantive, Bicap et Gampro sont des marques de commerce de Vantive Health Inc. ou de ses filiales.